

Illumina DRAGEN™ Bio-IT Platform

Raul Torrieri


Field Bioinformatics Scientist – Illumina LATAM

Setting the Stage



Investing in Secondary Analysis

Business Considerations for Investing in Secondary Analysis



01

Upfront Investments


Key things to consider when implementing secondary analysis include:

- Upfront hardware investments
- Bioinformatics & IT teams
- Compute and storage costs
- Direct instrument integration
- Security and compliance
- Ensuring accuracy
- Turnaround time
- Out of the box / Plug and play

Ability to Scale as You Grow

When adopting NGS, it's essential to choose a secondary analysis solution that is built for future expectations of scale.

02



Sequencing at Scale

Conducting NGS at scale poses a unique set of challenges



High Performance Secondary Analysis

01

Accurate

Detect SNPs and INDELS with high sensitivity and specificity

02

Ultra-rapid

Reduce analysis to ~25 minutes per 30x whole genome

03

Cost-efficient & scalable

Scale as needed while keeping costs low



Translating Sequencing Data into Insights




Library prep

Sequencing

Secondary analysis

Tertiary analysis

DRAGEN Germline Example



BCL
Conversion



Mapping
+ Aligning



Duplicate
Marking



Variant
Calling




From bases to variants in under 30 minutes



DRAGEN™ is Hardware-Accelerated Secondary Analysis


Dynamic Read Analysis for GENomics

DRAGEN Software Pipelines




Flexible Data Analysis

On-premise



BaseSpace™ Sequence Hub



Versatile Suite of Secondary Analysis Pipelines

All DRAGEN Pipelines can run on a single DRAGEN Platform

Illumina DRAGEN Germline Pipeline


End-to-end (BCL → VCF) NGS analysis, including advanced error model calibration for increased accuracy, and repeat expansion detection and genotyping through Illumina Expansion Hunter.

Illumina DRAGEN Somatic Pipeline

Includes tumor-only and tumor-normal modes, designed for detecting somatic variants in tumor samples.

DRAGEN RNA Pipeline

Performs transcriptome analysis starting with splice junction discovery and alignment, followed by gene fusion detection.



DRAGEN Joint Genotyping Pipeline

Calls variants across the genomes of multiple individuals, and scales to large cohorts of samples at expedited speeds with uncompromising accuracy.

DRAGEN CNV Pipeline

Performs copy number variant (CNV) analysis for germline and somatic exomes and genomes. Various levels of filtering can be applied to mitigate false positives before emitting the final calls

DRAGEN Map + Align Pipeline

Capable of ultra-rapid mapping and aligning DNA and RNA for both exomes and genomes.

New and enhanced DRAGEN pipelines are released periodically. For a comprehensive list, visit www.Illumina.com/DRAGEN

The Values of DRAGEN

illumina®



Easy to Implement

From command line to push-button

Bioinformatics Expertise

Push-button

- Simple Graphical User Interface (GUI)
- Managed service
- HIPAA* and GDPR compliance
- Workgroup capabilities
- Easy data sharing

Single command launch

- Easy-to-learn Linux based Command Line Interface (CLI)
- Simple command line execution

Advanced command line

- Script back-to-back jobs
- Make different configuration files for different applications


BaseSpace
SEQUENCE HUB

On-premise

*HIPAA compatibility applies in the US only with BaseSpace Enterprise

Record-Breaking Analysis Speed

Guinness World Records®
Fastest Genetic Diagnosis



Sample > Answer
19.5 Hours

In 2018, **Rady Children's Institute for Genomic Medicine** set the Guinness World Records® for Fastest Genetic Diagnosis leveraging the Illumina DRAGEN Bio-IT Platform.

Guinness World Records®
Fastest Analysis of 1,000 Genomes

1000_{WHG}



2hr25min

In 2017, **Children's Hospital of Philadelphia (CHOP)** set the Guinness World Records® for Fastest Analysis of 1,000 Genomes using the Illumina DRAGEN Bio-IT Platform in the cloud.

Accelerated Speeds

Analyze more in less time



Efficiency ratio

Analyze up to 57 genomes/day
vs ~2.5/day with standard solutions




**S4 Flow Cell
Demultiplexing**



**FASTQ > VCF
WHG @ 30x coverage**


Accelerated Speeds

DRAGEN on-premise v. BWA+GATK



Accelerated Speeds

DRAGEN on BaseSpace Sequence Hub v. BWA+GATK



Accurate Data

Detects small variants, copy number variants and structural variants
with high analytical **sensitivity** and **specificity**

precision**FDA** 


DRAGEN identified all 50 hidden variants
and ranked 1st in the following categories*

Hidden Variants	Indel Precision	Indel Recall	Indel F-Score	SNP F-score	SNP Recall
50/50	1 st	1 st	1 st	1 st	1 st

*Amongst entries that identified all hidden variants


Accurate Data—Small Variant Calling SNV

DRAGEN
detects fewer false
positive and false
negative SNV calls



Accurate Data—Small Variant Calling Indels

DRAGEN
detects fewer false
positive and false
negative indels



DRAGEN Metrics

The DRAGEN Platform Produces a Robust Portfolio of Metrics



Library Prep QC



Analysis QC



Demultiplexing



Duplicate Reads



Raw Data Processing

Cost-Effective



1 DRAGEN Server can **replace up to 30** traditional compute instances

Reduce Hardware Investment


Less compute, storage, power and other associated costs.



Reduce Storage footprint by **50%** using compression (CRAM)

Smaller Storage Footprint

Seamless data compression under-the-hood with CRAM output



\$5/Genome* on BaseSpace Sequence Hub

Low Cost Cloud Analysis

Most DRAGEN pipelines run at \$5 per sample on BaseSpace Sequence Hub

*Approximate cost. Varies based on input sample used. FASTQ to SNVs.

DRAGEN in Practice

illumina®



Popular Application Areas




Rare Genetic Disease

Reduces turnaround times required for genomic analysis, when fast results can be a critical factor.



Oncology

Easy-to-implement, cost effective and accelerated secondary analysis of whole exomes and genomes for cancer research.



Agrigenomics

The DRAGEN Platform is well suited to handle large samples and repetitive elements. With a provided reference, the DRAGEN Platform can analyze animals and plants of varying genomic complexities.




Population Genomics

The DRAGEN Platform can rapidly analyze sequenced samples, and accelerate reanalysis as computational tools improve over time.

NA12878 Titration INDEL Recall and Precision




Data generated by SickKids June 2019



GATK4 Hard Filtering vs. DRAGEN ROC InDels



Data generated by SickKids June 2019




DRAGEN on BaseSpace™ Sequence Hub

illumina®

4

Combining the Power of DRAGEN with the Flexibility of the Cloud

A hand is shown reaching out from the left side of the frame towards a wireframe globe. The globe is composed of a network of white lines and dots, representing a global network or data flow. The background is a dark blue gradient with a subtle pattern of white dots and lines, suggesting a digital or cloud environment.

DRAGEN on BaseSpace Sequence Hub

DRAGEN on BaseSpace Sequence Hub

Accurate, rapid secondary analysis in an easy-to-use, cloud-based environment



Accurate, Fast Analysis



Simple Workflow



Low-cost, Scalable Platform




Secure, Compliant Environment

Available Pipelines

- ✓ DRAGEN Germline Pipeline
- ✓ DRAGEN Somatic Pipeline
- ✓ DRAGEN Enrichment Pipeline (Q4 19)
- ✓ DRAGEN RNA Pipeline
- ✓ DRAGEN Joint Genotyping Pipeline
- ✓ DRAGEN Methylation Pipeline
- ✓ DRAGEN Reference Builder

High Performance Analysis, Simple Workflow



BaseSpace Sequence Hub—Security in the Cloud

Independently audited



Available with Enterprise Account



DRAGEN

Pipeline Overviews



Mapping & Aligning—Germline Pipeline Example




Mapping Aligning




Position Sorting*





Duplicate Marking*




Variant Calling



Position Sorting—Germline Pipeline Example



Duplicate Marking—Germline Pipeline Example



Variant Calling—Germline Pipeline Example



Mapping Aligning




Position Sorting*




Duplicate Marking*



Variant Calling





DRAGEN Germline Pipeline



DRAGEN Germline Pipeline—Accuracy


Highly accurate variant calling at low coverage



The DRAGEN v3 Germline Pipeline demonstrates high analytical sensitivity with relatively low false positive counts at 50x coverage. Performance comparable to GATK 4.1 at 50x coverage is possible with DRAGEN v3 at coverages lower than 50x.


DRAGEN Somatic Pipeline

Tumor/Normal Mode




DRAGEN Somatic Pipeline

Tumor-only Mode




DRAGEN RNA Pipeline

Gene Fusion Detection & Quantification



DRAGEN Joint Genotyping



DRAGEN Automatic QC Metrics Reporting: Mapper

Removes the need to run downstream tools for QC

- ✓ Number of samples
- ✓ Reads Processed
- ✓ Total
- ✓ Biallelic
- ✓ Multiallelic
- ✓ SNPs
- ✓ INDELs
- ✓ MNPs
- ✓ SNP Transitions
- ✓ SNP Transversions
- ✓ Ti/Tv ratio

- ✓ Heterozygous
- ✓ Homozygous
- ✓ Het/Hom ratio
- ✓ In dbSNP
- ✓ Novel
- ✓ Total
- ✓ Biallelic
- ✓ Multiallelic
- ✓ SNPs
- ✓ INDELs
- ✓ MNPs

DRAGEN Automatic QC Metrics Reporting: Variant Caller

- Total input reads
 - Number of duplicate reads (marked not removed)
 - Number of unique reads
 - Reads with mate sequenced
 - Reads without mate sequenced
 - QC-failed reads
 - Mapped reads
 - Number of unique & mapped reads (excl. dups)
 - Unmapped reads
 - Singleton reads (itself mapped; mate unmapped)
 - Paired reads (itself & mate mapped)
 - Properly paired reads
 - Not properly paired reads (discordant)
 - Reads with MAPQ [40:inf)
 - Reads with MAPQ [30:40)
 - Reads with MAPQ [20:30)
 - Reads with MAPQ [10:20)
 - Total reads in RG
 - Supplementary (chimeric) alignments
 - Average sequenced coverage over genome
- Total alignments
 - Secondary alignments
 - Supplementary (chimeric) alignments
 - Estimated read length
 - Bases in reference genome
 - Bases in target bed [% of genome]
 - Average sequenced coverage over genome
 - Average alignment coverage over genome
 - PCT of genome with coverage [40x:inf)
 - PCT of genome with coverage [30x:40x)
 - PCT of genome with coverage [20x:30x)
 - PCT of genome with coverage [10x:20x)
 - PCT of genome with coverage [5x:10x)
 - PCT of genome with coverage [2x: 5x)
 - PCT of genome with coverage [1x: 2x)
 - PCT of genome with coverage [0x: 1x)
 - DRAGEN mapping rate [mil. reads/second]
 - Secondary alignments
 - Estimated read length
 - Insert length: mean
- Number of duplicate reads (marked)
 - Number of unique reads
 - Reads with mate sequenced
 - Reads without mate sequenced
 - QC-failed reads
 - Mapped reads
 - Number of unique & mapped reads (excl. dups)
 - Unmapped reads
 - Singleton reads (itself mapped; mate unmapped)
 - Paired reads (itself & mate mapped)
 - Properly paired reads
 - Not properly paired reads (discordant)
 - Reads with MAPQ [40:inf)
 - Reads with MAPQ [30:40)
 - Reads with MAPQ [20:30)
 - Reads with MAPQ [10:20)
 - Reads with MAPQ [0:10)
 - Total alignments
 - Secondary alignments
 - Insert length: standard deviation

Thank You

techsupport@illumina.com | rtorrieri@illumina.com

illumina®

© 2020 Illumina, Inc. All rights reserved.